



IVFtech Sterile DS double sided cabinet

IVFtech DS Sterile cabinet is a double-sided cabinet which allows the embryologists to work from both sides simultaneously – providing a space saving and sociable solution. The cabinet can have up to 4 heated working areas with stereo microscopes and 2 extra small heated areas in the middle all independently controlled. Ideal for both small and larger laboratories as a peninsular or island unit.

Options

Sizes are available in: 120 and 180cm Monitor outside Options table top:

- Standard heat 450 x 650mm
- Extra heat 450 x 200mm
- Witness preparation flush
- IVFtech heated glass Ø45mm
- Table top integrated incubator

Filtration

Class H-14 HEPA filters in accordance with EN 1822. Filter efficiency 99.999% against 0.3 μm particle size. 99.995% in MPPS.

HiMacs or Stainless steel:

The Hi-Macs® table tops provides an excellent alternative to Stainless Steel. The Hi-Macs® surface is stain resistant and easy clean. Hi-Macs® is scratch resistant and doesn't scratch as easy as stainless steel, furthermore the Hi-Macs® surface can be refurbished and repolished. Easy to see the writing on the dish when working on the Hi-Macs® surface.

Advantages

- Easy to use
- Easy to service
- Easy to clean
- Hi-Macs® table tops
- Slope front for ergonomic work
- More work areas in less space
- Effective work flow



Maintenance & Service

All maintenance and service of filters can be performed from one front of the cabinet which in turn reduces down time during service.

Products to be used with IVFtech Sterile DS double sided cabinet

- Warming blocks
- Pull-out shelf
- Accessories for humidification
- Mini incubator
- Loud speakers

IVFtech

Toppevadvej 34-38 DK - 3660 Stenløse

Phone + 45 39 40 25 65 Fax + 45 39 40 25 64 www.ivftech.dk E-mail: info@ivftech.dk